




| Ø    | Optical access | Aperture X | ØT   | H1   | H2   | H3   | Ht   | Htc  | Z      |
|------|----------------|------------|------|------|------|------|------|------|--------|
| 3.00 | Ø1.2mm         | 30°        | 2.44 | 0.28 | 0.56 | 1.11 | 1.95 |      | 37°    |
| 3.50 |                |            | 2.94 | 0.28 | 0.56 |      | 1.95 |      | 32     |
| 4.00 | Ø1.9mm         | 60°        | 3.34 | 0.33 | 0.48 | 0.69 | 1.50 |      | 19°    |
| 4.00 | Ø1.9mm         | 50°        | 3.34 | 0.33 | 0.48 | 1.23 | 1.95 |      | 29°    |
| 4.00 | Ø1.9mm         | 40°        | 3.34 | 0.33 | 0.68 | 1.49 | 2.50 |      | 37°    |
| 4.00 |                |            | 3.80 | 0.10 | 0.44 | 0.66 | 1.42 | 1.20 | 23.75° |

**CONFIDENTIALITY NOTICE:** This document is for the sole use of the intended recipients and may contain confidential information. Any unauthorized review, use, disclosure or distribution is prohibited.

MDAC

DATE 23/04/2014

ECH. 1:10

Romain CIURKO

Matériau

Tolérance générale: +/- 0.025mm +/-0.15°

DESIGNATION: standard diamond anvil

**BETSA®**  
All for very high pressure experiments  
Z.I. 6 bis Rue de la Commune de Paris  
77 370 Nangis - FRANCE  
Tel: +33 1 64 08 27 33  
Fax: +33 1 64 08 33 51  
E mail: contact@betsa.fr